


SOLCO. The Trustworthy choice in your bag.

4CIS[®] PINEHURST

Anterior Cervical Plate System




OVERVIEW

The Pinehurst Anterior Cervical Plate System consists of plates and screws.

All components are manufactured with titanium alloy (Ti-6Al-4V).

The plates are available in lengths ranging from 08-91mm (Hole to hole) to accommodate on to five-level segments of fixation. The system also includes 3.5mm and 4.0mm titanium alloy screws which are available in lengths ranging from 10mm - 20mm

FEATURES & BENEFITS




DESIGN

Thickness : 2.0mm

Width

Screw Holes : 16.5mm

Neck : 14mm




Radius Curvature : Sagittal Plane


1 Level : 120mm

2 - 5 Levels : 200mm

Axial Plane : 25mm




Sagittal Plane




Axial Plane


Plate contouring




Locking CAP mechanism



Cephalad/caudal angulation



Medial/lateral angulation





Variable angle screws

- Cephalad/caudal: 23° range
- Medial/lateral: 10° range


Fixed angle screws

- Cephalad/caudal: Offset of 10°
- Medial/lateral: Offset of 10°



Screws

- Screws are color coded to identify diameter
- Regular screw diameter 3.5mm
- Each screw type is also available with diameter 4.0 mm for revision or where higher purchase is required



ONE LEVEL PLATES

Part Number	Level	Size(Hole-to-Hole)	Size(End-to-End)
8013-0108	One level plate	08mm	17mm
8013-0110	One level plate	10mm	19mm
8013-0112	One level plate	12mm	21mm
8013-0114	One level plate	14mm	23mm
8013-0116	One level plate	16mm	25mm
8013-0118	One level plate	18mm	27mm
8013-0120	One level plate	20mm	29mm
8013-0122	One level plate	22mm	31mm
8013-0124	One level plate	24mm	33mm
8013-0126	One level plate	26mm	35mm



TWO LEVEL PLATES

Part Number	Level	Size(Hole-to-Hole)	Size(End-to-End)
8013-0222	Two level plate	22mm	31mm
8013-0224	Two level plate	24mm	33mm
8013-0226	Two level plate	26mm	35mm
8013-0228	Two level plate	28mm	37mm
8013-0230	Two level plate	30mm	39mm
8013-0232	Two level plate	32mm	41mm
8013-0234	Two level plate	34mm	43mm
8013-0236	Two level plate	36mm	45mm
8013-0238	Two level plate	38mm	47mm
8013-0240	Two level plate	40mm	49mm
8013-0242	Two level plate	42mm	51mm
8013-0244	Two level plate	44mm	53mm
8013-0246	Two level plate	46mm	55mm




THREE LEVEL PLATES

Part Number	Level	Size(Hole-to-Hole)	Size(End-to-End)
8013-0336	Three level plate	36mm	45mm
8013-0339	Three level plate	39mm	48mm
8013-0342	Three level plate	42mm	51mm
8013-0345	Three level plate	45mm	54mm
8013-0348	Three level plate	48mm	57mm
8013-0351	Three level plate	51mm	60mm
8013-0354	Three level plate	54mm	63mm
8013-0357	Three level plate	57mm	66mm
8013-0360	Three level plate	60mm	69mm
8013-0363	Three level plate	63mm	72mm
8013-0366	Three level plate	66mm	75mm
8013-0369	Three level plate	69mm	78mm



FOUR, FIVE LEVEL PLATES

Part Number	Level	Size(Hole-to-Hole)	Size(End-to-End)
8013-0446	Four level plate	46mm	55mm
8013-0450	Four level plate	50mm	59mm
8013-0454	Four level plate	54mm	63mm
8013-0458	Four level plate	58mm	67mm
8013-0462	Four level plate	62mm	71mm
8013-0466	Four level plate	66mm	75mm
8013-0470	Four level plate	70mm	79mm
8013-0474	Four level plate	74mm	83mm
8013-0478	Four level plate	78mm	87mm
8013-0571	Five level plate	71mm	80mm
8013-0576	Five level plate	76mm	85mm
8013-0581	Five level plate	81mm	90mm
8013-0586	Five level plate	86mm	95mm
8013-0591	Five level plate	91mm	100mm



VARIABLE SELF DRILING SCREW

Part Number	Description	Size(Thread Length)
8208-3510	Self-Driling Ø3.5 Screw	10mm
8208-3511	Self-Driling Ø3.5 Screw	11mm
8208-3512	Self-Driling Ø3.5 Screw	12mm
8208-3513	Self-Driling Ø3.5 Screw	13mm
8208-3514	Self-Driling Ø3.5 Screw	14mm
8208-3515	Self-Driling Ø3.5 Screw	15mm
8208-3516	Self-Driling Ø3.5 Screw	16mm
8208-3517	Self-Driling Ø3.5 Screw	17mm
8208-3518	Self-Driling Ø3.5 Screw	18mm
8208-3519	Self-Driling Ø3.5 Screw	19mm
8208-3520	Self-Driling Ø3.5 Screw	20mm
8208-4010	Self-Driling Ø4.0 Screw	10mm
8208-4011	Self-Driling Ø4.0 Screw	11mm
8208-4012	Self-Driling Ø4.0 Screw	12mm
8208-4013	Self-Driling Ø4.0 Screw	13mm
8208-4014	Self-Driling Ø4.0 Screw	14mm
8208-4015	Self-Driling Ø4.0 Screw	15mm
8208-4016	Self-Driling Ø4.0 Screw	16mm
8208-4017	Self-Driling Ø4.0 Screw	17mm
8208-4018	Self-Driling Ø4.0 Screw	18mm
8208-4019	Self-Driling Ø4.0 Screw	19mm
8208-4020	Self-Driling Ø4.0 Screw	20mm

FIXED SELF DRILING SCREW

Part Number	Description	Size(Thread Length)
8206-3510	Self-Driling Ø3.5 Screw	10mm
8206-3511	Self-Driling Ø3.5 Screw	11mm
8206-3512	Self-Driling Ø3.5 Screw	12mm
8206-3513	Self-Driling Ø3.5 Screw	13mm
8206-3514	Self-Driling Ø3.5 Screw	14mm
8206-3515	Self-Driling Ø3.5 Screw	15mm
8206-3516	Self-Driling Ø3.5 Screw	16mm
8206-3517	Self-Driling Ø3.5 Screw	17mm
8206-3518	Self-Driling Ø3.5 Screw	18mm
8206-3519	Self-Driling Ø3.5 Screw	19mm
8206-3520	Self-Driling Ø3.5 Screw	20mm
8206-4010	Self-Driling Ø4.0 Screw	10mm
8206-4011	Self-Driling Ø4.0 Screw	11mm
8206-4012	Self-Driling Ø4.0 Screw	12mm
8206-4013	Self-Driling Ø4.0 Screw	13mm
8206-4014	Self-Driling Ø4.0 Screw	14mm
8206-4015	Self-Driling Ø4.0 Screw	15mm
8206-4016	Self-Driling Ø4.0 Screw	16mm
8206-4017	Self-Driling Ø4.0 Screw	17mm
8206-4018	Self-Driling Ø4.0 Screw	18mm
8206-4019	Self-Driling Ø4.0 Screw	19mm
8206-4020	Self-Driling Ø4.0 Screw	20mm

VARIABLE SELF TAPPING SCREW

Part Number	Description	Size(Thread Length)
8209-3510	Self-Tapping Ø3.5 Screw	10mm
8209-3511	Self-Tapping Ø3.5 Screw	11mm
8209-3512	Self-Tapping Ø3.5 Screw	12mm
8209-3513	Self-Tapping Ø3.5 Screw	13mm
8209-3514	Self-Tapping Ø3.5 Screw	14mm
8209-3515	Self-Tapping Ø3.5 Screw	15mm
8209-3516	Self-Tapping Ø3.5 Screw	16mm
8209-3517	Self-Tapping Ø3.5 Screw	17mm
8209-3518	Self-Tapping Ø3.5 Screw	18mm
8209-3519	Self-Tapping Ø3.5 Screw	19mm
8209-3520	Self-Tapping Ø3.5 Screw	20mm
8209-4010	Self-Tapping Ø4.0 Screw	10mm
8209-4011	Self-Tapping Ø4.0 Screw	11mm
8209-4012	Self-Tapping Ø4.0 Screw	12mm
8209-4013	Self-Tapping Ø4.0 Screw	13mm
8209-4014	Self-Tapping Ø4.0 Screw	14mm
8209-4015	Self-Tapping Ø4.0 Screw	15mm
8209-4016	Self-Tapping Ø4.0 Screw	16mm
8209-4017	Self-Tapping Ø4.0 Screw	17mm
8209-4018	Self-Tapping Ø4.0 Screw	18mm
8209-4019	Self-Tapping Ø4.0 Screw	19mm
8209-4020	Self-Tapping Ø4.0 Screw	20mm

FIXED SELF TAPPING SCREW

Part Number	Description	Size(Thread Length)
8207-3510	Self-Tapping Ø3.5 Screw	10mm
8207-3511	Self-Tapping Ø3.5 Screw	11mm
8207-3512	Self-Tapping Ø3.5 Screw	12mm
8207-3513	Self-Tapping Ø3.5 Screw	13mm
8207-3514	Self-Tapping Ø3.5 Screw	14mm
8207-3515	Self-Tapping Ø3.5 Screw	15mm
8207-3516	Self-Tapping Ø3.5 Screw	16mm
8207-3517	Self-Tapping Ø3.5 Screw	17mm
8207-3518	Self-Tapping Ø3.5 Screw	18mm
8207-3519	Self-Tapping Ø3.5 Screw	19mm
8207-3520	Self-Tapping Ø3.5 Screw	20mm
8207-4010	Self-Tapping Ø4.0 Screw	10mm
8207-4011	Self-Tapping Ø4.0 Screw	11mm
8207-4012	Self-Tapping Ø4.0 Screw	12mm
8207-4013	Self-Tapping Ø4.0 Screw	13mm
8207-4014	Self-Tapping Ø4.0 Screw	14mm
8207-4015	Self-Tapping Ø4.0 Screw	15mm
8207-4016	Self-Tapping Ø4.0 Screw	16mm
8207-4017	Self-Tapping Ø4.0 Screw	17mm
8207-4018	Self-Tapping Ø4.0 Screw	18mm
8207-4019	Self-Tapping Ø4.0 Screw	19mm
8207-4020	Self-Tapping Ø4.0 Screw	20mm

ORDERING INFORMATION INSTRUMENTS

Part Number	Description	Part Number	Description	Part Number	Description
5932-1117	Temporary Pin	5932-1124	AO I-Handle	9901-8004	Temporary Pin Caddy
5932-1118	Temporary Pin Insertor	5932-1121	Drill Guide-Fixed	9901-8007	ACP 1-3Level Plate Caddy
5932-1119	Self Constrained Awl	5932-1122	Drill Guide-Variable	9901-8003	ACP Screw Caddy
5932-1106	Drill 10mm	5932-1114	Plate Holder	9901-8008	ACP 4Level Plate Caddy
5932-1107	Drill 12mm	5932-1115	Plate Bender	9901-8009	ACP 5Level Plate Caddy
5932-1108	Drill 14mm	5932-1123	Screw Driver	9901-8012	Case
5932-1109	Drill 16mm	5932-1116	Hexlobe screw Removal tool		
5932-1211	Drill 18mm	5932-1111	Double DTS		
5932-1113	Tap 10mm	5932-1112	Single DTS		

4CIS[®] PINEHURST

SOLCO.

Solco Biomedical Corp.

5072 West Plano Parkway, Suite 210, Plano, TX 75093

Tel. 972-247-2486

Fax. 972-247-2413

E-mail: info@fg-solco.com

Website: www.solcospine.com

