

SOLCO. The Trustworthy choice in your bag.

4CIS[®] CHIRON


Spinal Fixation System




OVERVIEW

The Chiron Spinal Fixation System consists of a variety of pedicle screws, rods, nuts, cross links, and rod connectors. Implant components can be rigidly locked into a variety of different configurations to suit the individual pathology and anatomical conditions of the mature patient. All components are made of titanium alloy, per ASTM F136. From the thoracic spine to the ilium, the Chiron Spinal System facilitates surgeon choice and flexibility across patient types with a variety of implant options for treating multiple spinal pathologies with one system.

FEATURES & BENEFITS



SPECIFICATION



ORDERING INFORMATION

CHIRON STRAIGHT RODS

Cat.No	Description	Size
4302-5531	Spinal Rigid Rod	Ø5.5 30mm
4302-5536	Spinal Rigid Rod	35mm
4302-5504	Spinal Rigid Rod	40mm
4302-5501	Spinal Rigid Rod	45mm
4302-5505	Spinal Rigid Rod	50mm
4302-5555	Spinal Rigid Rod	55mm
4302-5506	Spinal Rigid Rod	60mm
4302-5565	Spinal Rigid Rod	65mm
4302-5507	Spinal Rigid Rod	70mm
4302-5575	Spinal Rigid Rod	75mm
4302-5508	Spinal Rigid Rod	80mm
4302-5585	Spinal Rigid Rod	85mm
4302-5509	Spinal Rigid Rod	90mm
4302-5595	Spinal Rigid Rod	95mm
4302-5510	Spinal Rigid Rod	100mm
4302-5511	Spinal Rigid Rod	110mm
4302-5512	Spinal Rigid Rod	120mm
4302-5513	Spinal Rigid Rod	130mm
4302-5514	Spinal Rigid Rod	140mm
4302-5515	Spinal Rigid Rod	150mm
4302-5516	Spinal Rigid Rod	160mm
4302-5517	Spinal Rigid Rod	170mm
4302-5518	Spinal Rigid Rod	180mm
4302-5519	Spinal Rigid Rod	190mm
4302-5520	Spinal Rigid Rod	200mm
4302-5525	Spinal Rigid Rod	250mm
4302-5530	Spinal Rigid Rod	300mm
4302-5540	Spinal Rigid Rod	400mm
4302-5550	Spinal Rigid Rod	500mm

CHIRON CROSSLINK

Cat.No	Description	Size
4412-3034	Crosslink	30~34mm
4412-3442	Crosslink	34~42mm
4412-4053	Crosslink	40~53mm
4412-5070	Crosslink	50~70mm

CHIRON PRE-BEND RODS

Cat.No	Description	Size
4332-5531	Spinal Pre-bend Rod	Ø5.5 30mm
4332-5536	Spinal Pre-bend Rod	35mm
4332-5504	Spinal Pre-bend Rod	40mm
4332-5545	Spinal Pre-bend Rod	45mm
4332-5505	Spinal Pre-bend Rod	50mm
4332-5555	Spinal Pre-bend Rod	55mm
4332-5506	Spinal Pre-bend Rod	60mm
4332-5565	Spinal Pre-bend Rod	65mm
4332-5507	Spinal Pre-bend Rod	70mm
4332-5575	Spinal Pre-bend Rod	75mm
4332-5508	Spinal Pre-bend Rod	80mm
4332-5509	Spinal Pre-bend Rod	90mm
4332-5510	Spinal Pre-bend Rod	100mm
4332-5511	Spinal Pre-bend Rod	110mm
4332-5512	Spinal Pre-bend Rod	120mm
4332-5513	Spinal Pre-bend Rod	130mm
4332-5514	Spinal Pre-bend Rod	140mm
4332-5515	Spinal Pre-bend Rod	150mm
4332-5516	Spinal Pre-bend Rod	160mm
4332-5517	Spinal Pre-bend Rod	170mm
4332-5518	Spinal Pre-bend Rod	180mm
4332-5519	Spinal Pre-bend Rod	190mm
4332-5520	Spinal Pre-bend Rod	200mm

CHIRON NUT SCREW

Cat.No	Description	Size
4205-0001	Chiron Nut Screw	M10 4.8mm

CHIRON ROD CONNECTOR

Cat.No	Description	Size
4502-0011	Chiron Rod Connector	Axial
4502-1001	Chiron Rod Connector	Domino

CHIRON ILIAC ROD CONNECTOR

Cat.No	Description	Size
4862-5515	Opened Connector	15mm
4862-5520	Opened Connector	20mm
4862-5525	Opened Connector	25mm
4862-5530	Opened Connector	30mm
4862-5540	Opened Connector	40mm
4862-5550	Opened Connector	50mm
4862-5560	Opened Connector	60mm
4863-5515	Opened Connector Variable	15mm
4863-5520	Opened Connector Variable	20mm
4863-5525	Opened Connector Variable	25mm
4863-5530	Opened Connector Variable	30mm
4863-5540	Opened Connector Variable	40mm
4863-5550	Opened Connector Variable	50mm
4863-5560	Opened Connector Variable	60mm
4864-5515	Closed Connector	15mm
4864-5520	Closed Connector	20mm
4864-5525	Closed Connector	25mm
4864-5530	Closed Connector	30mm
4864-5540	Closed Connector	40mm
4864-5550	Closed Connector	50mm
4864-5560	Closed Connector	60mm
4865-5515	Closed Connector Variable	15mm
4865-5520	Closed Connector Variable	20mm
4865-5525	Closed Connector Variable	25mm
4865-5530	Closed Connector Variable	30mm
4865-5540	Closed Connector Variable	40mm
4865-5550	Closed Connector Variable	50mm
4865-5560	Closed Connector Variable	60mm
4866-5515	Closed Connector Variable Curved	15mm
4866-5520	Closed Connector Variable Curved	20mm
4866-5525	Closed Connector Variable Curved	25mm
4866-5530	Closed Connector Variable Curved	30mm
4866-5540	Closed Connector Variable Curved	40mm
4866-5550	Closed Connector Variable Curved	50mm
4866-5560	Closed Connector Variable Curved	60mm

ORDERING INFORMATION INSTRUMENTS

Item.No	Description	Qty	Remark	Item.No	Description	Qty	Remark
4901-8028	4CIS® Awl	1		4901-8103	4CIS® I Ratcheting Handle	2	
9807-0038	Cannulated Awl		Option	4901-8104	4CIS® T Ratcheting Handle	2	
4901-0074	Guide Wire		Option	4901-8051	4CIS® French Rod Bender	1	
29-10012	Pedicle Probe for 5.5mm Screw	1		29-10026	Rod Pusher	1	
4901-0041	Pedicle Probe for 6.5mm Screw	1		4901-8012	4CIS® Chiron Rod Rocker 5.5 type	1	
4901-8001	4CIS® Pedicle Probe, Straight	1		4901-8058	4CIS® Bone Screw Driver I Handle	1	
4901-8002	4CIS® Pedicle Probe, Curved	1		4901-0055	Compressor	1	
29-10015	Guide Pins, Triangle	6		4901-0056	Distractor	1	
29-10016	Guide Pins, Ellipse	6		9807-0052	Torque Limit Handle (Nut)	1	
4901-8003	4CIS® Sounder, Straight	1		4901-7003	Nut Guide Cutting Forceps	1	
4901-8004	4CIS® Sounder, Angled	1		4901-0030	Screw Driver 3HEX	1	
4901-8005	4CIS® Chiron Tap Ø4.5	1		29-10022	Rod Gripper	1	
4901-8006	4CIS® Chiron Tap Ø5.5	1		4901-8093	4CIS® Chiron Persuader Tower type		Option
4901-8007	4CIS® Chiron Tap Ø6.5	1		4901-8109	4CIS® Chiron Persuader I type (Forceps)	1	
4901-8008	4CIS® Chiron Tap Ø7.5	1		4901-8117	4CIS® Chiron Persuader I type (Sleeve)	1	
4901-8110	4CIS® Chiron Tap Ø4.5 Cannulated		Option	4901-8097	4CIS® Chiron Nut Starter T-Handle Long type		Option
4901-8073	4CIS® Chiron Tap Ø5.5 Cannulated		Option	4901-8099	4CIS® Chiron Nut Final Driver Starter type		Option
4901-8074	4CIS® Chiron Tap Ø6.5 Cannulated		Option	4901-8105	4CIS® Nut Screw Final Driver Shaft 240mm	1	
4901-8075	4CIS® Chiron Tap Ø7.5 Cannulated		Option	4901-8140	4CIS® Nut Screw Final Driver Shaft 240mm	1	
4901-8094	4CIS® Chiron Screw Driver	2		4901-8139	4CIS® Chiron Anti-torque Wrench	1	
4901-8135	4CIS® Chiron Reduction Screw Driver	2		9901-4043	Rod Caddy	1	
4901-8048	4CIS® Nut Starter	2		9901-4044	Accessory Caddy	1	
4901-8013	4CIS® Nut Starter T Handle	2		9901-4049	Nut Screw Caddy	1	
4901-8031	4CIS® Screw Head Positioner	1		9901-4040	Implant Case	1	
4901-8032	4CIS® Chiron Rod Holder A type	1		9901-4041	Instrument Case 1	1	
4901-8049	4CIS® Chiron Rod Holder B Type	1		9901-4042	Instrument Case 2	1	

4CIS[®]
CHIRON

SOLCO.

Solco Biomedical Corp.

5072 West Plano Parkway, Suite 210, Plano, TX 75093

Tel. 972-247-2486

Fax. 972-247-2413

E-mail: info@fg-solco.com

Website: www.solcospine.com

