

SOLCO. The Trustworthy choice in your bag.

4CIS[®]

CHIRON XT


Minimally Invasive Spine Surgery system




OVERVIEW

Minimally invasive procedure is performed using K-wires and fluoroscopy, which allows smaller incisions rather than one long midline incision. Straight-forward instrumentation, self-tapping screws, and break-away tulip designs help simplifying the procedure and reduce operative time.

FEATURES & BENEFITS



SPECIFICATION



ORDERING INFORMATION

POLY CANULATED SCREW

Part No.	Description	Size	
4150-4520	Poly Cannulated Screw	Ø4.5	20mm
4150-4525	Poly Cannulated Screw	Ø4.5	25mm
4150-4530	Poly Cannulated Screw	Ø4.5	30mm
4150-4535	Poly Cannulated Screw	Ø4.5	35mm
4150-4540	Poly Cannulated Screw	Ø4.5	40mm
4150-4545	Poly Cannulated Screw	Ø4.5	45mm
4150-4550	Poly Cannulated Screw	Ø4.5	50mm
4150-5525	Poly Cannulated Screw	Ø5.5	25mm
4150-5530	Poly Cannulated Screw	Ø5.5	30mm
4150-5535	Poly Cannulated Screw	Ø5.5	35mm
4150-5540	Poly Cannulated Screw	Ø5.5	40mm
4150-5545	Poly Cannulated Screw	Ø5.5	45mm
4150-5550	Poly Cannulated Screw	Ø5.5	50mm
4150-5555	Poly Cannulated Screw	Ø5.5	55mm
4150-5560	Poly Cannulated Screw	Ø5.5	60mm
4150-6530	Poly Cannulated Screw	Ø6.5	30mm
4150-6535	Poly Cannulated Screw	Ø6.5	35mm
4150-6540	Poly Cannulated Screw	Ø6.5	40mm
4150-6545	Poly Cannulated Screw	Ø6.5	45mm
4150-6550	Poly Cannulated Screw	Ø6.5	50mm
4150-6555	Poly Cannulated Screw	Ø6.5	55mm
4150-6560	Poly Cannulated Screw	Ø6.5	60mm
4150-6565	Poly Cannulated Screw	Ø6.5	65mm
4150-6570	Poly Cannulated Screw	Ø6.5	70mm
4150-6575	Poly Cannulated Screw	Ø6.5	75mm
4150-6580	Poly Cannulated Screw	Ø6.5	80mm
4150-7530	Poly Cannulated Screw	Ø7.5	30mm
4150-7535	Poly Cannulated Screw	Ø7.5	35mm
4150-7540	Poly Cannulated Screw	Ø7.5	40mm
4150-7545	Poly Cannulated Screw	Ø7.5	45mm
4150-7550	Poly Cannulated Screw	Ø7.5	50mm
4150-7555	Poly Cannulated Screw	Ø7.5	55mm
4150-7560	Poly Cannulated Screw	Ø7.5	60mm
4150-7565	Poly Cannulated Screw	Ø7.5	65mm
4150-7570	Poly Cannulated Screw	Ø7.5	70mm
4150-7575	Poly Cannulated Screw	Ø7.5	75mm
4150-7580	Poly Cannulated Screw	Ø7.5	80mm
4150-8535	Poly Cannulated Screw	Ø8.5	35mm
4150-8540	Poly Cannulated Screw	Ø8.5	40mm
4150-8545	Poly Cannulated Screw	Ø8.5	45mm
4150-8550	Poly Cannulated Screw	Ø8.5	50mm
4150-8555	Poly Cannulated Screw	Ø8.5	55mm
4150-8560	Poly Cannulated Screw	Ø8.5	60mm
4150-8565	Poly Cannulated Screw	Ø8.5	65mm
4150-8570	Poly Cannulated Screw	Ø8.5	70mm
4150-8575	Poly Cannulated Screw	Ø8.5	75mm
4150-8580	Poly Cannulated Screw	Ø8.5	80mm

POLY CANULATED SCREW(HALF CLOSED TYPE)

Part No.	Description	Size	
4151-4520	Poly Cannulated Screw(Half Closed type)	Ø4.5	20mm
4151-4525	Poly Cannulated Screw(Half Closed type)	Ø4.5	25mm
4151-4530	Poly Cannulated Screw(Half Closed type)	Ø4.5	30mm
4151-4535	Poly Cannulated Screw(Half Closed type)	Ø4.5	35mm
4151-4540	Poly Cannulated Screw(Half Closed type)	Ø4.5	40mm
4151-4545	Poly Cannulated Screw(Half Closed type)	Ø4.5	45mm
4151-4550	Poly Cannulated Screw(Half Closed type)	Ø4.5	50mm
4151-5525	Poly Cannulated Screw(Half Closed type)	Ø5.5	25mm
4151-5530	Poly Cannulated Screw(Half Closed type)	Ø5.5	30mm
4151-5535	Poly Cannulated Screw(Half Closed type)	Ø5.5	35mm
4151-5540	Poly Cannulated Screw(Half Closed type)	Ø5.5	40mm
4151-5545	Poly Cannulated Screw(Half Closed type)	Ø5.5	45mm
4151-5550	Poly Cannulated Screw(Half Closed type)	Ø5.5	50mm
4151-5555	Poly Cannulated Screw(Half Closed type)	Ø5.5	55mm
4151-5560	Poly Cannulated Screw(Half Closed type)	Ø5.5	60mm
4151-6530	Poly Cannulated Screw(Half Closed type)	Ø6.5	30mm
4151-6535	Poly Cannulated Screw(Half Closed type)	Ø6.5	35mm
4151-6540	Poly Cannulated Screw(Half Closed type)	Ø6.5	40mm
4151-6545	Poly Cannulated Screw(Half Closed type)	Ø6.5	45mm
4151-6550	Poly Cannulated Screw(Half Closed type)	Ø6.5	50mm
4151-6555	Poly Cannulated Screw(Half Closed type)	Ø6.5	55mm
4151-6560	Poly Cannulated Screw(Half Closed type)	Ø6.5	60mm
4151-6565	Poly Cannulated Screw(Half Closed type)	Ø6.5	65mm
4151-6570	Poly Cannulated Screw(Half Closed type)	Ø6.5	70mm
4151-6575	Poly Cannulated Screw(Half Closed type)	Ø6.5	75mm
4151-6580	Poly Cannulated Screw(Half Closed type)	Ø6.5	80mm
4151-7530	Poly Cannulated Screw(Half Closed type)	Ø7.5	30mm
4151-7535	Poly Cannulated Screw(Half Closed type)	Ø7.5	35mm
4151-7540	Poly Cannulated Screw(Half Closed type)	Ø7.5	40mm
4151-7545	Poly Cannulated Screw(Half Closed type)	Ø7.5	45mm
4151-7550	Poly Cannulated Screw(Half Closed type)	Ø7.5	50mm
4151-7555	Poly Cannulated Screw(Half Closed type)	Ø7.5	55mm
4151-7560	Poly Cannulated Screw(Half Closed type)	Ø7.5	60mm
4151-7565	Poly Cannulated Screw(Half Closed type)	Ø7.5	65mm
4151-7570	Poly Cannulated Screw(Half Closed type)	Ø7.5	70mm
4151-7575	Poly Cannulated Screw(Half Closed type)	Ø7.5	75mm
4151-7580	Poly Cannulated Screw(Half Closed type)	Ø7.5	80mm
4151-8535	Poly Cannulated Screw(Half Closed type)	Ø8.5	35mm
4151-8540	Poly Cannulated Screw(Half Closed type)	Ø8.5	40mm
4151-8545	Poly Cannulated Screw(Half Closed type)	Ø8.5	45mm
4151-8550	Poly Cannulated Screw(Half Closed type)	Ø8.5	50mm
4151-8555	Poly Cannulated Screw(Half Closed type)	Ø8.5	55mm
4151-8560	Poly Cannulated Screw(Half Closed type)	Ø8.5	60mm
4151-8565	Poly Cannulated Screw(Half Closed type)	Ø8.5	65mm
4151-8570	Poly Cannulated Screw(Half Closed type)	Ø8.5	70mm
4151-8575	Poly Cannulated Screw(Half Closed type)	Ø8.5	75mm
4151-8580	Poly Cannulated Screw(Half Closed type)	Ø8.5	80mm

ORDERING INFORMATION

PRE-BEND ROD

Part No.	Description	Size	
4343-5515	Sharp Tip Pre-bend Rod	Ø5.5	15mm
4343-5520	Sharp Tip Pre-bend Rod	Ø5.5	20mm
4343-5525	Sharp Tip Pre-bend Rod	Ø5.5	25mm
4343-5530	Sharp Tip Pre-bend Rod	Ø5.5	30mm
4343-5535	Sharp Tip Pre-bend Rod	Ø5.5	35mm
4343-5540	Sharp Tip Pre-bend Rod	Ø5.5	40mm
4343-5545	Sharp Tip Pre-bend Rod	Ø5.5	45mm
4343-5550	Sharp Tip Pre-bend Rod	Ø5.5	50mm
4343-5555	Sharp Tip Pre-bend Rod	Ø5.5	55mm
4343-5560	Sharp Tip Pre-bend Rod	Ø5.5	60mm
4343-5565	Sharp Tip Pre-bend Rod	Ø5.5	65mm
4343-5570	Sharp Tip Pre-bend Rod	Ø5.5	70mm
4343-5575	Sharp Tip Pre-bend Rod	Ø5.5	75mm
4343-5580	Sharp Tip Pre-bend Rod	Ø5.5	80mm
4343-5585	Sharp Tip Pre-bend Rod	Ø5.5	85mm
4343-5590	Sharp Tip Pre-bend Rod	Ø5.5	90mm
4343-5595	Sharp Tip Pre-bend Rod	Ø5.5	95mm
4343-5100	Sharp Tip Pre-bend Rod	Ø5.5	100mm
4343-5110	Sharp Tip Pre-bend Rod	Ø5.5	110mm
4343-5120	Sharp Tip Pre-bend Rod	Ø5.5	120mm
4343-5130	Sharp Tip Pre-bend Rod	Ø5.5	130mm
4343-5140	Sharp Tip Pre-bend Rod	Ø5.5	140mm
4343-5150	Sharp Tip Pre-bend Rod	Ø5.5	150mm
4343-5160	Sharp Tip Pre-bend Rod	Ø5.5	160mm
4343-5170	Sharp Tip Pre-bend Rod	Ø5.5	170mm
4343-5180	Sharp Tip Pre-bend Rod	Ø5.5	180mm
4343-5190	Sharp Tip Pre-bend Rod	Ø5.5	190mm
4343-5200	Sharp Tip Pre-bend Rod	Ø5.5	200mm
4343-5250	Sharp Tip Pre-bend Rod	Ø5.5	250mm

STRAIGHT ROD

Part No.	Description	Size	
4344-5515	Sharp Tip Straight Rod	Ø5.5	15mm
4344-5520	Sharp Tip Straight Rod	Ø5.5	20mm
4344-5525	Sharp Tip Straight Rod	Ø5.5	25mm
4344-5530	Sharp Tip Straight Rod	Ø5.5	30mm
4344-5535	Sharp Tip Straight Rod	Ø5.5	35mm
4344-5540	Sharp Tip Straight Rod	Ø5.5	40mm
4344-5545	Sharp Tip Straight Rod	Ø5.5	45mm
4344-5550	Sharp Tip Straight Rod	Ø5.5	50mm
4344-5555	Sharp Tip Straight Rod	Ø5.5	55mm
4344-5560	Sharp Tip Straight Rod	Ø5.5	60mm
4344-5565	Sharp Tip Straight Rod	Ø5.5	65mm
4344-5570	Sharp Tip Straight Rod	Ø5.5	70mm
4344-5575	Sharp Tip Straight Rod	Ø5.5	75mm
4344-5580	Sharp Tip Straight Rod	Ø5.5	80mm
4344-5585	Sharp Tip Straight Rod	Ø5.5	85mm
4344-5590	Sharp Tip Straight Rod	Ø5.5	90mm
4344-5595	Sharp Tip Straight Rod	Ø5.5	95mm
4344-5100	Sharp Tip Straight Rod	Ø5.5	100mm
4344-5110	Sharp Tip Straight Rod	Ø5.5	110mm
4344-5120	Sharp Tip Straight Rod	Ø5.5	120mm
4344-5130	Sharp Tip Straight Rod	Ø5.5	130mm
4344-5140	Sharp Tip Straight Rod	Ø5.5	140mm
4344-5150	Sharp Tip Straight Rod	Ø5.5	150mm
4344-5160	Sharp Tip Straight Rod	Ø5.5	160mm
4344-5170	Sharp Tip Straight Rod	Ø5.5	170mm
4344-5180	Sharp Tip Straight Rod	Ø5.5	180mm
4344-5190	Sharp Tip Straight Rod	Ø5.5	190mm
4344-5200	Sharp Tip Straight Rod	Ø5.5	200mm
4344-5250	Sharp Tip Straight Rod	Ø5.5	250mm
4344-5300	Sharp Tip Straight Rod	Ø5.5	300mm
4344-5350	Sharp Tip Straight Rod	Ø5.5	350mm
4344-5400	Sharp Tip Straight Rod	Ø5.5	400mm

CHIRON NUT SCREW

Part No.	Description	Size	
4205-0001	Chiron Nut Screw	M10	4.8mm

ORDERING INFORMATION INSTRUMENTS

INSTRUMENTS

Item No.	Description	Qty	Remark
9807-0053	VP Needle Diamond type Long	4	
4901-3019	VP Needle Diamond type short		Option
9807-0038	4CIS® Cannulated Awl	1	
4901-0074	4CIS® Guide Wire	10	480mm
4901-8033	4CIS® Guide Wire Tool	1	
4901-8034	4CIS® Dilator 1	2	
4901-8035	4CIS® Dilator 2	2	
4901-8036	4CIS® Dilator 3	2	
4901-8037	4CIS® Chiron Tap 4.5 mm(Cannulated)		Option
4901-8038	4CIS® Chiron Tap 5.5 mm(Cannulated)	1	
4901-8039	4CIS® Chiron Tap 6.5 mm(Cannulated)	1	
4901-8040	4CIS® Chiron Tap 7.5 mm(Cannulated)	1	
4901-8027	4CIS® Inline Rasp		Option
4901-8041	4CIS® Ring C type	6	
4901-8042	4CIS® Closed Ring		Option
4901-8095	4CIS® Chiron XT Screw Driver	2	
4901-8103	4CIS® I Ratcheting Handle	2	
4901-8104	4CIS® T Ratcheting Handle	2	
4901-8107	4CIS® Chiron-XT Rod Length Gauge	1	
4901-8051	4CIS® French Rod Bender	1	
4901-8024	4CIS® Chiron-XT Rigid Rod Inserter	2	
4901-8046	4CIS® Rod Holder Forceps Type	1	
4901-8031	4CIS® Screw Head Positioner	1	Rod Checker
4901-8059	4CIS® Chiron XT Nut Starter	2	
4901-8060	4CIS® Chiron-XT Nut Starter T-handle	2	
4901-8097	4CIS® Nut Starter T-handle Long type		Option
4901-8025	4CIS® Chiron-XT Persuader	2	
4901-8065	4CIS® Chiron-XT Persuader (Half Closed Head)	2	For Half Closed Type
4901-8058	4CIS® Poly Axial Bone Screw Driver I Handle	1	
4901-8079	4CIS® Chiron XT Compressor&Distractor	1	
4901-8092	4CIS® Chiron XT Compression&Distraction Bar		Option
4901-8080	4CIS® Chiron XT Fixing Sleeve 1	2	
4901-8081	4CIS® Chiron XT Fixing Sleeve 2	2	
4901-8044	4CIS® Torque Stabilizer	1	
4901-8082	4CIS® Chiron XT Counter Torque	1	For Half Closed Type
4901-8108	4CIS® Chiron XT Counter Torque	1	For Open Type
4901-8061	4CIS® Chiron XT Nut Screw Final Driver	2	
9807-0052	Torque Limit Handle	1	
4901-8101	4CIS® Chiron XT Tab Breaker	1	
4901-8111	4CIS® Ring Cutter	1	For Half Closed Type
9901-4045	Rod Caddy	1	
9901-4047	Implant Case	1	
9901-4048	Instrument Case	1	

4CIS[®] CHIRON XT

SOLCO.

Solco Biomedical Corp.

5072 West Plano Parkway, Suite 210, Plano, TX 75093

Tel. 972-247-2486

Fax. 972-247-2413

E-mail: info@fg-solco.com

Website: www.solcospine.com

