


4CIS[®]

PEEK Interbody Fusion Cage

PEEK Interbody Fusion Cage

PEEK has the better modulus of elasticity to the cancellous bone than any other materials, such as the metal or the carbon fiber



A Wide Fusion Area


Comparing with other cages, it has a wide fusion area which guarantees the high union rate.



A serrated design to prevent cage sliding



Easy to couple and uncouple



A Lordotic configuration (0°, 4°, 8°)



Easy insertion



Fits the Anatomy



The bullet type cage has an anatomical similarity with the endplate.








PEEK Interbody Fusion Cage Implants List


Cat. No	Description	Size
 5800-0008	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×8mm, 0°
5800-0009	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×9mm, 0°
5800-0010	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×10mm, 0°
 5800-0011	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×11mm, 0°
5800-0012	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×12mm, 0°
5800-0013	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×13mm, 0°
5800-0014	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×14mm, 0°

Cat. No	Description	Size
 5804-0009	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×9mm, 4°
5804-0010	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×10mm, 4°
5804-0011	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×11mm, 4°
 5804-0012	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×12mm, 4°
5804-0013	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×13mm, 4°
5804-0014	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×14mm, 4°

Cat. No	Description	Size
 5808-0009	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×9mm, 8°
5808-0010	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×10mm, 8°
5808-0011	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×11mm, 8°
 5808-0012	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×12mm, 8°
5808-0013	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×13mm, 8°
5808-0014	PEEK Posterior Lumbar Interbody Fusion Cage	25×11×14mm, 8°

PEEK Interbody Fusion Cage Instruments List

Cat No.	Description	Size
 5901-0035	Reamer Sharp	0° × 7mm
5901-1038	Reamer Sharp	0° × 8mm
5901-0036	Reamer Sharp	0° × 9mm
5901-1039	Reamer Sharp	0° × 10mm
5901-0037	Reamer Sharp	0° × 11mm
5901-1040	Reamer Sharp	0° × 12mm
5901-0038	Reamer Sharp	0° × 13mm
5901-1041	Reamer Sharp	0° × 14mm

Cat No.	Description	Size
 5901-1033	Reamer Sharp	4° × 7mm
5901-1042	Reamer Sharp	4° × 8mm
5901-1034	Reamer Sharp	4° × 9mm
5901-1043	Reamer Sharp	4° × 10mm
5901-1035	Reamer Sharp	4° × 11mm
5901-1044	Reamer Sharp	4° × 12mm
5901-1036	Reamer Sharp	4° × 13mm
5901-1045	Reamer Sharp	4° × 14mm

* The design of products can be amended without prior consent.

PEEK Interbody Fusion Cage Instruments List

Cat No.	Description	Size
5901-0040	Reamer Blunt	7mm
5901-1046	Reamer Blunt	8mm
5901-0041	Reamer Blunt	9mm
5901-1047	Reamer Blunt	10mm
5901-0042	Reamer Blunt	11mm
5901-1048	Reamer Blunt	12mm
5901-0043	Reamer Blunt	13mm
5901-1049	Reamer Blunt	14mm

Cat No.	Description	Size
5901-0032	Bone Curette	Angled / Left
5901-0033	Bone Curette	Angled / Right

Cat No.	Description	Size
5901-0029	Osteotome	Angled / Left
5901-0030	Osteotome	Angled / Right

Cat No.	Description
5901-0045	T-Handle

Cat No.	Description
5901-0050	PLIF Last Impactor

Cat No.	Description
5901-0026	PEEK PLIF Cage Holder

Cat No.	Description
5901-0047	Graft Impactor

Cat No.	Description	Size
29-17014	Duramer Retractor	16×215mm

Cat No.	Description
70-17004	4CIS® Cage System Container

Cat No.	Description
5901-0046	Graft Holder

* The design of products can be amended without prior consent.